

# Recent Advances in Potentiometric Scanning Electrochemical Microscopy

author:  
András Kiss

supervisor:  
prof. Géza Nagy


Department of General and Physical Chemistry  
University of Pécs, Hungary



April 18, 2017



# Potentiometric Scanning Electrochemical Microscopy

## A Scanning Probe Microscopic technique



# Ion-selective micropipettes

As SECM probes





Valinomycin

$$E = E^\theta + \frac{RT}{z_i F} \ln \left[ a_i + \sum_j \left( k_{ij} a_j^{z_i/z_j} \right) \right]$$

Nikolsky-equation


# The problem with potentiometric SECM

Distortion at high scan rate



# Why is the image distorted?

The RC time constant



The time required to discharge the capacitor by  $\approx 37\%$  ( $1/e$ ).

Or to charge it by  $\approx 63\%$  ( $1 - 1/e$ ).

$$\tau = R \cdot C$$


$$R = 5 \text{ G}\Omega$$

$$C = 500 \text{ pF}$$

$$\tau = 2.5 \text{ s}$$


# Distortion of potentiometric imaging

In the case of a linescan




# Why is it so important to complete the scan quickly?

Example: corrosion of a magnesium alloy



Corrosion of the AZ63 magnesium-aluminium-zinc alloy.

## Trade-off triangle of potentiometric SECM




Solution #1: Solid contact micropipettes as SECM probes.


# Liquid vs. solid contact micropipettes

## Comparison of construction

Liquid contact




Solid contact



# Comparison of the electrodes' resistance


## Experimental setup and result



Type	$R_{ISME} / G\Omega$
Liquid contact	4.80
Solid contact	0.56



# Comparison of the electrodes' performance

## Experimental setup




# Comparison of the electrodes' performance

## Results



# Application in corrosion science: galvanic corrosion of Mg


## Experimental setup




# Application in corrosion science: galvanic corrosion of Mg

## Results

Liquid contact







Solid contact



Solution #2: Optimizing scanning patterns and algorithms.



# "Traditional" scanning algorithms for the 2D raster pattern

## 1. "Meander" and "fast comb" scanning algorithms




# "Traditional" scanning algorithms for the 2D raster pattern

## 2. "Comb" scanning algorithm




# The RC time constant




# New SECM scanning patterns based on the polar-coordinate system

Cartesian coordinate  
system based scanning  
pattern







Polar coordinate  
system based scanning  
pattern



# New SECM scanning patterns based on the polar-coordinate system

Simulated scans with the "web" and the "arc" algorithms



## Comparison of the simulated scans




Parameters of the simulations:

- 2000  $\mu\text{m} \times 2000 \mu\text{m}$  scan area,
- 100  $\mu\text{m} \times 100 \mu\text{m}$  resolution,
- 1 s for each data acquisition point,
- 500  $\mu\text{m}/\text{s}$  probe positioning speed.

Algorithm	#n	time (s)	mean squared error
Meander	441	440	$2.75 \times 10^{-2}$
Fast comb	441	520	$2.07 \times 10^{-2}$
Comb	441	881	$2.75 \times 10^{-2}$
Web	110	109	$9.63 \times 10^{-3}$
Arc	341	340	$2.95 \times 10^{-3}$

# Confirmation with experimental SECM scans


## Traditional scanning algorithms



Images recorded with the meander, fast comb, and comb algorithms.

# Confirmation with experimental SECM scans



## New scanning algorithms




Images recorded with the web, and the arc algorithms.

Scans are completed 75% and 20% faster, images have almost 10 times less distortion.

# Example 1: SECM study on the effectiveness of corrosion protection coatings using the arc scanning pattern




## Example 2: SECM study on the galvanic corrosion of the AZ63 alloy and iron, using the arc scanning pattern



pH-dependent potential map above the corroding (A) magnesium, and (B) iron samples.


Solution #2: Signal processing.

## Deconvolution of the distorted image



$$E_{cell}(t) = E_{cell}(\infty) + [E_{cell}(0) - E_{cell}(\infty)]e^{-t/RC}$$


$$E_{cell}(\infty) = \frac{[E_{cell}(t) - E_{cell}(0)]e^{-t/RC}}{1 - e^{-t/RC}}$$



Response of an antimony/antimony-oxid - Ag/AgCl/3M KCl cell when pH changes from 4 to 6. RC time constant was determined from the fit.

# Minimal working example:

## Deconvolution of a response to activity step





Transient response of the antimony microelectrode to analyte activity step (red), and deconvolutions performed with different  $RC$  time-constants. The measured time-constant was  $\tau = 3.76\text{s}$  (cyan).

## Comparison of the deconvoluted time-potential recordings.

$e^{-0.5/RC}$	$RC(s)$	Mean squared error
raw recording (0)	raw recording (0)	53.43
0.7	1.4	22.03
0.75	1.74	15.88
0.8	2.24	9.01
<b>0.8755</b>	<b>3.76</b>	<b>3.83</b>
0.9	4.75	16.99
0.95	9.75	781.94

# One step further



## Deconvolution of linescans




(A) Raw, and, (B) deconvoluted SECM linescans above the center of the target, at  $h = 100$  mm height, using three different equilibration periods: 0.5 s, (green), 1 s, (blue), and 5 s, (red).

# Deconvolution of the distorted image


pH microscopy with antimony microelectrode



# Linescans with magnesium-ion selective micropipettes, and their deconvolution




# 2D scans with magnesium-ion selective micropipettes, and their deconvolution




# Example 3: SECM study on the galvanic corrosion of a zinc - copper galvanic pair

With antimony microelectrode



## Example 4: Raw, and deconvoluted SECM image and microphoto of a corroding carbon-steel sample.

With antimony microelectrode



# Summary

- Compromise between image quality and scan time in potentiometric SECM.
- I've demonstrated two approaches to alleviate this compromise.
- With optimized scanning patterns and algorithms, scans are faster up to 4 times, imaging distortion is 10 times less.
- These algorithms work only with radially symmetric targets.
- Another approach is to deconvolution. Can be applied to any system. Image quality is drastically improved, and scan time is shorter.

**Thank you for your  
kind attention!**



**Pécs, Hungary**